

ERR-2	HL70357 ERR-3	HL70516 ERR-4	HL70533 ERR-5	ERR # Error Description ERR-5 ERR-5	Error Details ERR-8
	207	E	3	00 Duplicate message detected	Message had been ignored by ICARE system. Duplicate message detected
	207	E	4	00 Message had been Ignored	Message had been ignored by ICARE system (ICARE Error)
	207	E	4	01 Patient record is too small	Patient record is too small
	207	E	4	03 Invalid ID Type Code	Invalid ID Type Code
PID^3^X	207	E	4	04 Invalid ID	Invalid ID # X
RXA^X^3	207	E	5	10 Invalid CVX or NDC code	Invalid CVX or NDC code # X
OBX^X^Y^5	207	E	5	11 Invalid Adverse Event	Invalid Adverse Event Shot # X Oservation # Y
OBX^X^Y^5	207	E	5	12 Invalid Contraindication	Invalid Contraindication Shot # X Oservation # Y
OBX^X^Y^5	207	E	5	13 Invalid Immunity	Invalid Immunity Shot # X Oservation # Y
	207	E	4	14 Alias Name is empty	Alias Name is empty. Alias Record Ignored.
	207	E	4	15 Invalid Alias Code	Invalid Alias Code. Alias Record Ignored.
RXA^X	207	E	5	16 Invalid Origin Code	Invalid Origin Code for Shot # X
RXA^X^3	207	E	2	17 Invalid Shot Date	Invalid Shot # X Date
RXA^X^5	207	E	2	18 Invalid Shot Code	Invalid Shot Code # X
RXA^X	207	E	4	19 Invalid Shot Record	Invalid Shot Record # X
OBX^X^Y^5	207	E	5	20 Invalid Immunity	Invalid Immunity Shot # X Oservation # Y
PID^5	207	E	4	21 Invalid Patient Name	Invalid Patient Name
RXA^X^3	207	E	2	22 Invalid Shot Date	Invalid Shot # X Date
	101	E	4	0074 Wrong user name and/or password	Message had been Ignored. Wrong user name and/or password
QPD	101	E	4	0098 Invalid QPD Message	Message had been Ignored. Invalid QPD Message.
	207	E	4	0105 Internal Parser Error	Message had been Ignored.
MSH^4^1	101	E	4	0135 Empty Clinic Identifier	Message had been Ignored. Empty Clinic Identifier
MSH^4^1	101	E	4	0167 Invalid HL7 Clinic id and/or user id	HL7 Clinic id and/or user id are not present in ICARE TEST system
MSH^4^1	101	E	4	0178 Invalid HL7 Clinic id and/or user id	HL7 Clinic id and/or user id are not present in ICARE DEV system
MSH^4^1	101	E	4	0206 Invalid HL7 Clinic id and/or user id	HL7 Clinic id and/or user id are not present in ICARE
QRD	101	E	4	0245 Missing QRD Segment	QRD Segment is required
	207	E	4	0272 Internal Database Error	Internal Database Error. Please call customer support
QRD^STA	101	E	4	0333 QRD Segment Invalidate	QRD Segment must include filter STA (check last downloaded updates date)
RCP^2	101	E	4	0381 Too many records to select	Message had been Ignored. You can't select more than 35 records at once
VXQ^2	101	E	4	0579 Too many records to select	Message had been Ignored. You can't select more than 35 records at once
VXQ^QRD	101	E	4	0587 QRD Segment is missing	QRD Segment is required
VXQ^QRF	101	E	4	0595 QRF Segment is missing	QRF Segment is required
MSH^9^3	200	E	4	0734 Unsupported Message Type	X is not supported
	207	E	4	0741 Unspecified application internal error	Unspecified application internal error TEXT
VXU	101	E	4	0754 Message Too Large	X HL7 messages were submitted. Please submit one and only one HL7 message per transaction
MSH	101	E	4	0760 Empty message	Empty or invalid message"
	207	E	4	0774 Internal Parser Error	Message had been Ignored
RXA^X^5	101	E	4	0887 Invalid CVX/NDC code	Patient Shot Record # X - Shot CVX and NDC Codes are Empty or Invalid. Shot record had been Ignored.
RXA^X^3	101	E	4	0893 Invalid Administered Date	Patient Shot Record # X - Shot Date is Empty, Invalid or in the Future. Shot record had been Ignored.
	207	E	4	997 Parser Error	(Parser Error)
	207	E	4	1052 Internal DMBS Error	Internal DBMS Error (DBMS Error)
PID	207	E	4	1104 Load demographic error	Message had been Ignored. Load demographic error - ParserError
PID^7	101	E	2	1118 Invalid Birth Date	Message had been Ignored. Invalid Birth Date
PID^3^1	101	E	4	1124 Invalid Medical Record Number	Message had been Ignored. Medical Record Number is invalid
PID^5	101	E	4	1130 Invalid Patient Name	Message had been Ignored. Patient Name (First and/or Last) are Empty

PID	207	E	4	1162 Load demographic error	Message had been Ignored. Load demographic error - ParserError
PID^3^X	207	E	4	1206 Invalid Patient Id #	Patient Id Record # X had been Ignored. (ParserError)
PID^19	207	E	4	1226 Invalid Social Security #	Social Security # in PID.19 is invalid. Error - ParserError
PID^11^X	207	E	4	1273 Invalid Patient Address	Patient Address Record # X had been Ignored. Error - ParserError
PID^13^X	207	E	4	1301 Invalid Patient Telecom	Patient Telecom Record # X - Value is Empty or Invalid. Record had been Ignored
PID^13^X	207	E	4	1312 Invalid Patient Telecom	Patient Home Phone/Email Record # X had been Ignored. Error - ParserError
NK1^X^2	101	E	4	1384 Invalid Next of Kin Name	Next of Kin Record # X had been Ignored. Patient Name (First and/or Last) is Empty.
NK1^X^3	101	E	4	1384 Invalid Next of Kin Relationship code	Next of Kin Record # X had been Ignored. Patient Relationship Code is Empty
NK1^X^4	101	E	4	1458 Invalid Next of Kin Address	Next of Kin Record # X Address # 1 had been Ignored. Error - ParserError
NK1^X^Y^5	207	E	4	1494 Invalid Next of Kin Phone/E-mail	Telecom Record # Y of Next of Kin Record # X had been Ignored. Error - ParserError
NK1^X^Y^6	207	E	4	1527 Invalid NK Business Phone/E-mail	Business Telecom Record # Y of Next of Kin Record # X had been Ignored. Error - ParserError
NK1^X	207	E	4	1535 Invalid Next of Kin Segment	Next of Kin Record # X had been Ignored. Error - ParserError
	207	E	4	2160 Internal Database Error	Internal Database Error. Please call customer support
	0	I		01 Message Accepted	Message had been accepted by ICARE system
	207	I		999 Parser Information	(Parser Information)
PID^8	207	W		02 Sex code is undefined	Sex code is undefined
	207	W		05 Invalid Address Type Code	Invalid Address Type Code
PID^13^X	207	W		06 Invalid Telecom Location Code	Invalid Telecom # X Location Code
PID^13^X	207	W		07 Invalid Telecom Type Code	Invalid Telecom # X Type Code
NK1^X^3	207	W		08 Invalid NK1 Relationship Code	Invalid NK1 Relationship Code
NK1^X^4	207	W		09 Invalid Address Type Code	Invalid Address # X Type Code
RXA^X^15	207	W		10 I-CARE performed FUZZY search on inventory	I-CARE performed FUZZY search on inventory shot # X
RXA	207	W		0930 Duplicate Shot Segment	Patient Shot Record # X - Duplicate with previous record. Shot record had been Ignored.
OBX^5	101	W		0957 Invalid Event Code	Patient Shot Record # X with Event Record # Y - Event Code and Description are Empty.
OBX^3	101	W		0965 Invalid Observation ID	Patient Shot Record # X with Event Record # Y - Observation ID is Empty. Event Record had been Ignored.
OBX^3	201	W		0971 Event code not supported	Patient Shot Record # X with Event Record # Y - Observation ID is not supported
OBX^X^14	101	W		0978 Invalid event date	Patient Shot Record # X with Event Record # Y - Invalid Empty or in the Future Event Date
OBX^X^Y	101	W		998 Invalid OBX Segment	Shot Record # X with Adverse Event Record # Y had been Ignored. Error - ParserError
OBX^X	101	W		1006 Invalid OBX Segment	Shot Record # X with Adverse Event Record # Y had been Ignored. Error - GenericParserError
PID^3^X^5	101	W	4	1177 Invalid MRN Type	MRN # X Type is Empty
PID^3^X^1	101	W	4	1184 Invalid MRN	MRN # X is Empty or Invalid
PID^14^X	207	W	4	1339 Invalid Patient Telecom	Patient Telecom Record # X - Value is Empty or Invalid. Record had been Ignored
PID^14^X	207	W	4	1350 Invalid Patient Telecom	Patient Home Phone/Email Record # X had been Ignored. Error - (ParserError)